

FILM APPLICATION CHART

<p>Cell / Tissue / Bacterial Culture Fluorescence / Storage HPLC / Chromatography Luminescence / Microscopy Protein Crystallography Sonication</p>		Plate - Tube		Application						Product Attributes				Material									
		Raised-Rim Micro-Plates	Flat-Top Micro-Plates	PCR, 8-Tube Strips/Single Plate Row	Fluorescence / Storage	Luminescence / Storage	ELISA/Incubation/General	Cell/Tissue/Bact. Culture	Long Term Storage	Automation Compatible	Protein Crystallography	HPLC/Chromatography	Sonication	RNAse & DNase Free	Chemically Resistant	Pierceable	Optically Clear	Transparent	Light Blocking	Breathable	Temp. Range	Film/Foil Thickness	Adhesive Thickness
202537	AbsorbMax™ Black Sealing Film, Non-Sterile	✓		✓																	-40°C to +80°C	Vinyl 203µm / 8 mil	Acrylic 37µm / 1.45mil
202504	AeraSeal™ Sealing Film for Culture Plates, Non-Sterile	✓					✓							✓							-20°C to +80°C	Rayon / Acrylic	139.7 µm / 5.5mil
202505	AeraSeal™ Sealing Film for Culture Plates, Sterile	✓					✓							✓							-20°C to +80°C	Rayon / Acrylic	139.7 µm / 5.5mil
202504Roll	AeraSeal™ Sealing Film, Roll, Non-Sterile	✓					✓							✓							-20°C to +80°C	Rayon / Acrylic	139.7 µm / 5.5mil
202550NS	AeraSeal™ SealMate™ Dispenser Kit, Non-Sterile	✓					✓							✓							-20°C to +80°C	Rayon / Acrylic	139.7 µm / 5.5mil
202550	AeraSeal™ SealMate™ Dispenser Kit, Sterile	✓					✓							✓							-20°C to +80°C	Rayon / Acrylic	139.7 µm / 5.5mil
202551NS	AeraSeal™ Rolls for SealMate™ Dispenser, Non-Sterile	✓					✓							✓							-20°C to +80°C	Rayon / Acrylic	139.7 µm / 5.5mil
202551	AeraSeal™ Rolls for SealMate™ Dispenser, Sterile	✓					✓							✓							-20°C to +80°C	Rayon / Acrylic	139.7 µm / 5.5mil
202512	BrightMax™ Film, Non-Sterile	✓					✓											✓			-40°C to +80°C	Vinyl 203µm / 8mil	Acrylic 37µm / 1.45mil
202547	ThermalSeal RTS™ for RT qPCR, Optically Clear, Non-Sterile	✓	✓						✓	✓	✓			✓	✓	✓					-80°C to +105°C	Polyolefin 50.8µm / 2mil	Silicone 50.8µm / 2mil
202547S	ThermalSeal RTS™ for RT qPCR, Optically Clear, Sterile	✓	✓						✓	✓	✓			✓	✓	✓					-80°C to +105°C	Polyolefin 50.8µm / 2mil	Silicone 50.8µm / 2mil
202540	X-Pierce™ Film, Non-Sterile	✓							✓	✓	✓			✓		✓					-40°C to +90°C	Vinyl 88.9µm / 3.5mil	Acrylic 17.8µm / 0.7mil
202541	X-Pierce™ Film, Sterile	✓							✓	✓	✓			✓		✓					-40°C to +90°C	Vinyl 88.9µm / 3.5mil	Acrylic 17.8µm / 0.7mil
202515	Zone-Free™ Adhesive Film, Non-Sterile	✓							✓	✓				✓		✓					-40°C to +90°C	Polyethylene 69.9µm / 2.75mil	Acrylic 26.7µm / 1.05mil

Trademark - Excel Scientific, Inc.

FILM APPLICATION CHART

ELISA / Incubation / General Purpose

No.	Description	Plate - Tube		Application							Product Attributes					Material						
		Raised-Rim Micro-Plates	Flat-Top Micro-Plates	PCR & Tube Strips/ Single Plate Row / Strip Wells	qPCR (Real-Time)	PCR	ELISA/Incubation/General Cell/Tissue/Bact. Culture	Long Term Storage	Automation Compatible	Protein Crystallography	HPLC/Chromatography	Sonication	Temp. Range	RNAse & DNAse Free	Chemically Resistant	Pierceable	Optically Clear	Transparent	Light Blocking	Breathable	Film/Foil Thickness	Adhesive Thickness
202557	EZcap™ General Purpose FilmStrips™, Non-Sterile	✓	✓			✓							-40°C to +120°C	✓			✓				Polypropylene 50.8µm / 2mil	Acrylic 24.1µm / 0.95mil
202558	EZcap™ General Purpose FilmStrips™, Sterile	✓	✓			✓							-40°C to +120°C	✓			✓				Polypropylene 50.8µm / 2mil	Acrylic 24.1µm / 0.95mil
202527	SealPlate® ColorTab™ Film, Non-Sterile, Assorted Colors (Blue, Green, Lavender, Orange, Red), Non-Sterile	✓				✓							-40°C to +120°C				✓				Polyester 50.81µm / 2mil	Acrylic 26.7µm / 1.05mil
202502	SealPlate® Sealing Film, Clear, Non-Sterile	✓				✓							-40°C to +120°C				✓				Polyester 50.8µm / 2mil	Acrylic 26.7µm / 1.05mil
202503	SealPlate® Sealing Film, Clear, Sterile	✓				✓							-40°C to +120°C				✓				Polyester 50.8µm / 2mil	Acrylic 26.7µm / 1.05mil
202552NS	SealPlate® SealMate™ Dispenser Kit, Non-Sterile	✓				✓							-40°C to +120°C				✓				Polyester 50.8µm / 2mil	Acrylic 26.7µm / 1.05mil
202552	SealPlate® SealMate™ Dispenser Kit, Sterile	✓				✓							-40°C to +120°C				✓				Polyester 50.8µm / 2mil	Acrylic 26.7µm / 1.05mil
202553NS	SealPlate® Rolls for SealMate™ Dispenser Kit, Non-Sterile	✓				✓							-40°C to +120°C				✓				Polyester 50.8µm / 2mil	Acrylic 26.7µm / 1.05mil
202553	SealPlate® Rolls for SealMate™ Dispenser Kit, Sterile	✓				✓							-40°C to +120°C				✓				Polyester 50.8µm / 2mil	Acrylic 26.7µm / 1.05mil
202506	SealPlate® MiniStrips™, Non-Sterile	✓	✓			✓							-40°C to +120°C				✓				Polyester 50.8µm / 2mil	Acrylic 26.7µm / 1.05mil
202507	SealPlate® MiniStrips™, Sterile	✓	✓			✓							-40°C to +120°C				✓				Polyester 50.8µm / 2mil	Acrylic 26.7µm / 1.05mil
202522RM	SealPlate® RollMate™ Dispenser w/Pre-Scored Sheets, Non-Sterile	✓				✓							-40°C to +120°C				✓				Polyester 50.8µm / 2mil	Acrylic 26.7µm / 1.05mil
202511NS	ThinSeal™ Film, Non-Sterile	✓				✓							-40°C to +120°C				✓				Polyester 25.4µm / 1mil	Acrylic 24.1µm / 0.95mil
202511	ThinSeal™ Film, Sterile	✓				✓							-40°C to +120°C				✓				Polyester 25.4µm / 1mil	Acrylic 24.1µm / 0.95mil

Trademark - Excel Scientific, Inc.

FILM APPLICATION CHART

Automation Compatible
Long Term Storage

No.	Description	Plate - Tube		Application							Product Attributes					Material					
		Raised-Rim Micro-Plates	Flat-Top Micro-Plates	qPCR	PCR	Luminescence / Microscopy	Cell/Tissue/Bact. Culture	Long Term Storage	Automation Compatible	Protein Crystallography	HPLC/Chromatography	Sonication	Temp. Range	RNAse & DNase Free	Chemically Resistant	Pierceable	Optically Clear	Transparent	Light Blocking	Breathable	Film/Foil Thickness
202535	AlumaSeal CS™ Sealing Film, Non-Sterile	✓					✓					-80°C to +120°C	✓			✓				Aluminum Foil 50.8µm / 2mil	Acrylic 50.8µm / 2mil
202536	AlumaSeal CS™ Sealing Film, Sterile	✓					✓					-80°C to +120°C	✓			✓				Aluminum Foil 50.8µm / 2mil	Acrylic 50.8µm / 2mil
248897	AlumaSeal II™ Adhesive Aluminum Sealing Foil, Non-Sterile	✓			✓		✓	✓				-80°C to +120°C	✓	✓		✓				Aluminum Foil 35.6µm / 1.4mil	Acrylic 27.9µm / 1.1mil
248898	AlumaSeal II™ Adhesive Aluminum Sealing Foil, Sterile	✓			✓		✓	✓				-80°C to +120°C	✓	✓		✓				Aluminum Foil 35.6µm / 1.4mil	Acrylic 27.9µm / 1.1mil
248897Roll	AlumaSeal® Adhesive Aluminum Sealing Foil, Non-Sterile, Roll	✓			✓		✓	✓				-80°C to +120°C	✓	✓		✓				Aluminum Foil 38.1µm / 1.5mil	Acrylic 37.6µm / 1.48mil
202548	eXTReme™ FoilSeal™, Non-Sterile	✓	✓				✓	✓				-80°C to +120°C	✓	✓	✓	✓				Aluminum Foil 48.3µm / 1.9mil	Acrylic 30.5µm / 1.2mil
202549	eXTReme™ FoilSeal™, Large, Non-Sterile	✓	✓				✓	✓				-80°C to +120°C	✓	✓	✓	✓				Aluminum Foil 48.3µm / 1.9mil	Acrylic 30.5µm / 1.2mil
202519	EZ Pierce Adhesive Sealing Film, Non-Sterile	✓						✓				-40°C to +90°C	✓	✓						Polyethylene 69.9µm / 2.75mil	Acrylic 24.1µm / 0.95mil
202520	EZ Pierce Adhesive Sealing Film, Sterile	✓						✓				-40°C to +90°C	✓	✓						Polyethylene 69.9µm / 2.75mil	Acrylic 24.1µm / 0.95mil
202555	EZCap™ FoilStrips™, Non-Sterile	✓	✓		✓		✓	✓				-80°C to +120°C	✓	✓	✓	✓				Aluminum Foil 50.8µm / 2mil	Acrylic 36.8µm / 1.45mil
202510	FoilSeal™ RollMate™ Film Dispenser with Pre-Scored Sheets	✓	✓		✓		✓	✓				-80°C to +120°C	✓	✓	✓	✓				Aluminum Foil 50.8µm / 2mil	Acrylic 36.8µm / 1.45mil
202547	ThermalSeal RTS™ for RT qPCR, Optically Clear, Non-Sterile	✓	✓		✓		✓	✓	✓			-80°C to +105°C	✓	✓	✓					Polyolefin 50.8µm / 2mil	Silicone 50.8µm / 2mil
202547S	ThermalSeal RTS™ for RT qPCR, Optically Clear, Sterile	✓	✓		✓		✓	✓	✓			-80°C to +105°C	✓	✓	✓					Polyolefin 50.8µm / 2mil	Silicone 50.8µm / 2mil
252000Roll	ThermalSeal RTS™ for RT qPCR, Optically Clear, Non-Sterile, Roll	✓	✓		✓		✓	✓	✓			-80°C to +105°C	✓	✓	✓					Polyolefin 50.8µm / 2mil	Silicone 50.8µm / 2mil
202540	X-Pierce™ Film, Non-Sterile	✓					✓	✓	✓			-40°C to +90°C		✓	✓					Vinyl 88.9µm / 3.5mil	Acrylic 17.8µm / 0.7mil
202541	X-Pierce™ Film, Sterile	✓					✓	✓	✓			-40°C to +90°C		✓	✓					Vinyl 88.9µm / 3.5mil	Acrylic 17.8µm / 0.7mil
202515	Zone-Free™ Adhesive Film, Non-Sterile	✓					✓	✓				-40°C to +90°C		✓	✓					PP/PE 69.9µm / 2.75mil	Acrylic 26.7µm / 1.05mil

Trademark - Excel Scientific, Inc.

FILM APPLICATION CHART

PCR and **qPCR**

No.	Description	Plate - Tube			Application							Product Attributes					Material					
		Raised-Rim Micro-Plates	Flat-Top Micro-Plates	PCR & Tube Strips/Single Plate Row	qPCR (Real-Time)	PCR	ELISA/Incubation/General	Cell/Tissue/Bact. Culture	Long Term Storage	Automation Compatible	Protein Crystallography	HPLC/Chromatography	Sonication	Temp. Range	RNAse & DNase Free	Chemically Resistant	Pierceable	Optically Clear	Transparent	Light Blocking	Breathable	Film/Foil Thickness
248897	AlumaSeal II® Adhesive Aluminum Sealing Foil, Non-Sterile	✓			✓		✓	✓					-80°C to +120°C	✓	✓		✓				Aluminum Foil 35.6µm / 1.4mil	Acrylic 27.9µm / 1.1mil
248898	AlumaSeal II® Adhesive Aluminum Sealing Foil, Sterile	✓			✓		✓	✓					-80°C to +120°C	✓	✓		✓				Aluminum Foil 35.6µm / 1.4mil	Acrylic 27.9µm / 1.1mil
248897ROLL	AlumaSeal® Adhesive Aluminum Sealing Foil, Non-Sterile, Roll	✓			✓		✓	✓					-80°C to +120°C	✓	✓		✓				Aluminum Foil 37.8µm / 1.5mil	Acrylic 37.8µm / 1.48mil
202560	eXTreme™ Seal, Clear, Non-Sterile	✓	✓		✓	✓							-20°C to +120°C	✓		✓					Aluminum Foil 48.3µm / 1.9mil	Acrylic 30.5µm / 1.2mil
202561	eXTreme™ Seal, Large, Clear, Non-Sterile	✓			✓	✓							-20°C to +120°C	✓		✓					Polyester 50.8µm / 2mil	Acrylic 41.9µm / 1.65mil
202555	EZCap™ FoilStrips™, Non-Sterile	✓	✓		✓		✓	✓					-80°C to +120°C	✓	✓	✓		✓			Polyester 50.8µm / 2mil	Acrylic 41.9µm / 1.65mil
202556	EZCap™ PCR FilmStrips™, Non-Serile	✓	✓		✓	✓							-40°C to +120°C	✓		✓					Aluminum Foil 50.8µm / 2mil	Acrylic 36.8µm / 1.45mil
202510	FoilSeal™ RollMate™ Film Dispenser with Pre-Scored Sheets	✓	✓		✓		✓	✓					-80°C to +120°C	✓	✓	✓		✓			Polyester 50.8µm / 2mil	Acrylic 41.9µm / 1.65mil
202500	ThermalSeal® Film, Clear, Non-Sterile	✓			✓								-40°C to +120°C	✓	✓		✓				Aluminum Foil 50.8µm / 2mil	Acrylic 36.8µm / 1.45mil
202501	ThermalSeal® Film, Clear, Sterile	✓			✓								-40°C to +120°C	✓	✓		✓				Polypropylene 50.8µm / 2mil	Acrylic 24.1µm / 0.95mil
202509	ThermalSeal® MiniStrips™, Non-Sterile	✓			✓								-40°C to +120°C	✓			✓				Polypropylene 50.8µm / 2mil	Acrylic 24.1µm / 0.95mil
202545	ThermalSeal A® Sealing Film, Non-Sterile	✓	✓		✓								-40°C to +120°C	✓			✓				Polypropylene 50.8µm / 2mil	Acrylic 24.1µm / 0.95mil
202547	ThermalSeal RTS™ for RT qPCR, Optically Clear, Non-Sterile	✓	✓		✓		✓	✓	✓				-80°C to +105°C	✓	✓	✓					Polyester 50.8µm / 2mil	Acrylic 41.9µm / 1.65mil
202547S	ThermalSeal RTS™ for RT qPCR, Optically Clear, Sterile	✓	✓		✓		✓	✓	✓				-80°C to +105°C	✓	✓	✓					Polyolefin 50.8µm / 2mil	Silicone 50.8µm / 2mil
252000	ThermalSeal RT™ for RT qPCR, Optically Clear, Non-Sterile	✓			✓								-40°C to +120°C	✓		✓					Polyolefin 50.8µm / 2mil	Silicone 50.8µm / 2mil
252005	ThermalSeal RT2RR™ for RT qPCR, Optically Clear, Non-Sterile	✓	✓		✓								-40°C to +120°C	✓		✓					Polyester 50.8µm / 2mil	Acrylic 29.2µm / 1.15mil
252000ROLL	ThermalSeal RTS™ for RT qPCR, Optically Clear, Non-Sterile, Roll	✓	✓		✓		✓	✓	✓				-80°C to +105°C	✓	✓	✓					Polyester 50.8µm / 2mil	Acrylic 29.2µm / 1.15mil
																					Polyolefin 50.8µm / 2mil	Silicone 50.8µm / 2mil

Trademark - Excel Scientific, Inc.